

## Tau mouse Monoclonal Antibody Biotin Conjugated

Catalog No: #C92901Bio

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Tau mouse Monoclonal Antibody Biotin Conjugated
Host Species	mouse
Clonality	Monoclonal
Purification	affinity purified by Protein G
Applications	WB IHC ICC IF
Species Reactivity	Hu Ms Rat
Conjugates	Bio
Target Name	Tau
Other Names	Tau
Accession No.	Uniprot:P10636Gene ID:4137
Uniprot	P10636
GeneID	4137
Calculated MW	52/79
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

WB=1:50-200 IHC-P=1:50-200 IHC-F=1:50-200 ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only