

# ATG5/APG5L mouse Monoclonal Antibody PE Conjugated

Catalog No: #C92904PE

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	ATG5/APG5L mouse Monoclonal Antibody PE Conjugated
Host Species	mouse
Clonality	Monoclonal
Purification	affinity purified by Protein G
Applications	ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	full length
Conjugates	PE
Target Name	ATG5/APG5L
Other Names	ATG5/APG5L
Accession No.	Uniprot:Q9H1Y0Gene ID:9474
Uniprot	Q9H1Y0
GeneID	9474
Calculated MW	32
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only