

ATG7/APG7 mouse Monoclonal Antibody PE Conjugated



Catalog No: #C92905PE

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	ATG7/APG7 mouse Monoclonal Antibody PE Conjugated
Host Species	mouse
Clonality	Monoclonal
Purification	affinity purified by Protein G
Applications	ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	71-333/703
Conjugates	PE
Target Name	ATG7/APG7
Other Names	ATG7/APG7
Accession No.	Uniprot:O95352Gene ID:10533
Uniprot	O95352
GeneID	10533
Calculated MW	78
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only