

P63 mouse Monoclonal Antibody PE Conjugated

Catalog No: #C92907PE

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	P63 mouse Monoclonal Antibody PE Conjugated
Host Species	mouse
Clonality	Monoclonal
Purification	affinity purified by Protein G
Applications	ICC
Species Reactivity	Hu Ms Rat
Conjugates	PE
Target Name	P63
Other Names	P63
Accession No.	Uniprot:Q9H3D4Gene ID:8626
Uniprot	Q9H3D4
GeneID	8626
Calculated MW	77
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

ICC=1:50-200

Note: This product is for in vitro research use only