

Myoglobin mouse Monoclonal Antibody AP Conjugated

Catalog No: #C92908AP

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Myoglobin mouse Monoclonal Antibody AP Conjugated
Host Species	mouse
Clonality	Monoclonal
Purification	affinity purified by Protein G
Applications	WB IHC ICC
Species Reactivity	Hu Rat
Conjugates	AP
Target Name	Myoglobin
Other Names	Myoglobin
Accession No.	Uniprot:P02144Gene ID:4151
Uniprot	P02144
GeneID	4151
Calculated MW	17
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:10000-20000 IHC-P=1:50-200 IHC-F=1:50-200 ICC=1:50-200

Note: This product is for in vitro research use only