

# Mouse 8-iso Prostaglandin (8-iso-PG) ELISA Kit

Catalog No: #EK20186



Package Size: #EK20186-1 48T #EK20186-2 96T

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse 8-iso Prostaglandin (8-iso-PG) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Ms

## Application Details

Detect Range:15.6-1000pg/ml

Sensitivity:5.80pg/ml

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 &mu;L

Assay Time:1-4.5h

Detection wavelength:450 nm

## Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only