

Mouse Nicotinamide adenine dinucleotide(NADH) Elisa kit

Catalog No: #EK20189

Orders: order@signalwayantibody.com

Package Size: #EK20189-1 48T #EK20189-2 96T

Support: tech@signalwayantibody.com

Description

Product Name	Mouse Nicotinamide adenine dinucleotide(NADH) Elisa kit
Brief Description	Elisa kit
Applications	ELISA
Species Reactivity	Ms

Application Details

Detect Range:62.5-4000pg/mL

Sensitivity:18.2pg/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

References

Detection Method:sandwich

Note: This product is for in vitro research use only