

# General Gamma-Aminobutyric Acid (gABA) ELISA Kit

Catalog No: #EK20300



Package Size: #EK20300-1 48T #EK20300-2 96T

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	General Gamma-Aminobutyric Acid (gABA) ELISA Kit
Applications	ELISA
Species Reactivity	General

## Application Details

Detect Range:31.2-2000pg/mL

Sensitivity:10.9pg/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 &mu;L

Assay Time:1-4.5h

Detection wavelength:450 nm

## Product Description

Detection Method:competitive inhibition

Note: This product is for in vitro research use only