

Mouse total bilirubin (TB) ELISA Kit

Catalog No: #EK20307



Package Size: #EK20307-1 48T #EK20307-2 96T

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mouse total bilirubin (TB) ELISA Kit
Applications	ELISA
Species Reactivity	Mouse

Application Details

Detect Range:0.78-50ng/mL

Sensitivity:0.28ng/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μ L

Assay Time:1-4.5h

Detection wavelength:450 nm

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only