

# Rat Glycerol-3-phosphate dehydrogenase, mitochondrial (GPD2) ELISA Kit



Catalog No: #EK20309

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Package Size: #EK20309-1 48T #EK20309-2 96T

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |  |
|--------------------|--|
| Product Name       | Rat Glycerol-3-phosphate dehydrogenase, mitochondrial (GPD2) ELISA Kit |
| Brief Description  | ELISA Kit  |
| Applications       | ELISA  |
| Species Reactivity | Rat  |

## Application Details

Detect Range:31.2-2000pg/mL

Sensitivity:12.4pg/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 &mu;L

Assay Time:1-4.5h

Detection wavelength:450 nm

## Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only