

Human Glutamate receptor delta-2 subunit (GRID2) ELISA Kit



Catalog No: #EK20246

Orders: order@signalwayantibody.com

Package Size: #EK20246-1 48T #EK20246-2 96T

Support: tech@signalwayantibody.com

Description

Product Name	Human Glutamate receptor delta-2 subunit (GRID2) ELISA Kit
Applications	ELISA
Species Reactivity	Hu

Application Details

Detect Range:0.156-10ng/mL

Sensitivity:0.067ng/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only