

PAH Rabbit mAb

Catalog No: #52041

Package Size: #52041-1 50ul #52041-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	PAH Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal antibody
Clone No.	S06-3H4
Isotype	Rabbit IgG
Purification	Affinity Purified
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	A synthetic peptide of human PAH
Conjugates	Unconjugated
Modification	Unmodification
Other Names	Phe-4-monooxygenase
Accession No.	Swiss-Prot:P00439GenelD:5053
Uniprot	P00439
GenelD	5053
Calculated MW	Calculated MW: 52 kDa; Observed MW: 52 kDa
Concentration	0.3 mg/ml
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

WB: 1/1000-1/5000;

Images



Western blot detection of PAH in 3T3,C6,Rat Brain,K562 cell lysates using PAH Rabbit mAb(1:1000 diluted).Predicted band size:52kDa.Observed band size:52kDa.

Background

Swiss-Prot Acc.P00439. Catalyzes the hydroxylation of L-phenylalanine to L-tyrosine.

Note: This product is for in vitro research use only