

Neurofascin Rabbit mAb

Catalog No: #52044

Package Size: #52044-1 50ul #52044-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	Neurofascin Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal antibody
Clone No.	S07-0B3
Isotype	Rabbit IgG
Purification	Affinity Purified
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant protein of mouse Neurofascin
Conjugates	Unconjugated
Modification	Unmodification
Other Names	Neurofascin
Accession No.	Swiss-Prot:Q810U3GeneID:269116
Uniprot	Q810U3
GeneID	269116
Calculated MW	Calculated MW: 150 kDa; Observed MW: 186,155,75 kDa
Concentration	0.3 mg/ml
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

WB: 1/1000

Images



Western blot detection of Neurofascin in Mouse Brain,Rat Brain,COS7,HeLa cell lysates using Neurofascin Rabbit mAb(1:1000 diluted).Predicted band size:150kDa.Observed band size:186,155,75kDa.

Background

Cell adhesion, ankyrin-binding protein which may be involved in neurite extension, axonal guidance, synaptogenesis, myelination and neuron-glia cell interactions.

Note: This product is for in vitro research use only