

CD90 Rabbit mAb

Catalog No: #52065

Package Size: #52065-1 50ul #52065-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	CD90 Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal antibody
Clone No.	S09-6G7
Isotype	Rabbit IgG
Purification	Affinity Purified
Applications	WB
Species Reactivity	Human,Rat
Immunogen Description	A synthetic peptide of human CD90
Conjugates	Unconjugated
Modification	Unmodification
Other Names	Thy-1 membrane glycoprotein; Thy-1 antigen; CDw90; CD90; THY1;
Accession No.	Swiss-Prot:P04216GenelD:7070
Uniprot	P04216
GenelD	7070
Calculated MW	Calculated MW: 18 kDa; Observed MW: 22 kDa
Concentration	0.3 mg/ml
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

WB: 1/1000;

Images



Western blot detection of CD90/Thy1 in Rat Brain cell using CD90/Thy1 Rabbit mAb(1:1000 diluted). Predicted band size:18kDa.Observed band size:22kDa.

Background

Swiss-Prot Acc.P04216.May play a role in cell-cell or cell-ligand interactions during synaptogenesis and other events in the brain.

Note: This product is for in vitro research use only