

Annexin V Rabbit mAb

Catalog No: #52083

Package Size: #52083-1 50ul #52083-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	Annexin V Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal antibody
Clone No.	S08-3A4
Isotype	Rabbit IgG
Purification	Affinity Purified
Applications	WB IF
Species Reactivity	Human,Mouse,Rat
Immunogen Description	A synthetic peptide of human Annexin V
Conjugates	Unconjugated
Modification	Unmodification
Other Names	PP4;ANX5;ENX2;RPRGL3.
Accession No.	Swiss-Prot:P08758GenelD:308
Uniprot	P08758
GenelD	308
Calculated MW	Calculated MW: 36 kDa; Observed MW: 36 kDa
Concentration	0.3 mg/ml
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details


WB: 1/2000; ICC/IF: 1/100;

Images



Western blot detection of Annexin V in K562,C6,3T3,HeLa cell lysates using Annexin V Rabbit mAb(1:1000 diluted).Predicted band size:36kDa.Observed band size:36kDa.

Immunofluorescence of Annexin V(green) in Hela cells using Annexin V Rabbit mAb at dilution 1/200, and DAPI(blue)



Background

Swiss-Prot Acc.P08758.This protein is an anticoagulant protein that acts as an indirect inhibitor of the thromboplastin-specific complex, which is involved in the blood coagulation cascade.

Note: This product is for in vitro research use only