

AAMP Rabbit mAb

Catalog No: #52113

Package Size: #52113-1 50ul #52113-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	AAMP Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal antibody
Clone No.	S01-8F4
Isotype	Rabbit IgG
Purification	Affinity Purified
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	A synthetic peptide of human AAMP
Conjugates	Unconjugated
Modification	Unmodification
Accession No.	Swiss-Prot:Q13685GeneID:14
Uniprot	Q13685
GeneID	14
Calculated MW	Calculated MW: 47 kDa; Observed MW: 47 kDa
Concentration	0.3 mg/ml
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

WB: 1/1000;

Images



Western blot detection of AAMP in C6,HeLa cell lysates using AAMP Rabbit mAb(1:1000 diluted).Predicted band size:47kDa.Observed band size:47kDa.

Background

Swiss-Prot Acc.Q13685.Plays a role in angiogenesis and cell migration. In smooth muscle cell migration, may act through the RhoA pathway.

Note: This product is for in vitro research use only