

ACADS Rabbit mAb

Catalog No: #52116

Package Size: #52116-1 50ul #52116-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	ACADS Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal antibody
Clone No.	S08-4C2
Isotype	Rabbit IgG
Purification	Affinity Purified
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	A synthetic peptide of human ACADS
Conjugates	Unconjugated
Modification	Unmodification
Other Names	SCAD; ACAD3
Accession No.	Swiss-Prot:P16219GenelD:35
Uniprot	P16219
GenelD	35
Calculated MW	Calculated MW: 44 kDa; Observed MW: 44 kDa
Concentration	0.3 mg/ml
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

WB: 1/1000;

Images



Western blot detection of ACADS in Jurkat,C6,HeLa cell lysates using ACADS Rabbit mAb(1:1000 diluted).Predicted band size:44kDa.Observed band size:44kDa.

Background

Swiss-Prot Acc.P16219.Introduces a double bond at position 2 in saturated acyl-CoA's of short chain length, i.e. less than 6 carbon atoms.

Note: This product is for in vitro research use only