

ACAT1 Rabbit mAb

Catalog No: #52117

Package Size: #52117-1 50ul #52117-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	ACAT1 Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal antibody
Clone No.	S02-5G8
Isotype	Rabbit IgG
Purification	Affinity Purified
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	A synthetic peptide of human ACAT1
Conjugates	Unconjugated
Modification	Unmodification
Other Names	T2; MAT; ACAT; THIL
Accession No.	Swiss-Prot:P24752GenelD:38
Uniprot	P24752
GenelD	38
Calculated MW	Calculated MW: 45 kDa; Observed MW: 45 kDa
Concentration	0.3 mg/ml
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

WB: 1/1000;

Images



Western blot detection of ACAT1 in K562, HeLa cell lysates using ACAT1 Rabbit mAb (1:1000 diluted). Predicted band size: 45 kDa. Observed band size: 45 kDa.

Background

Swiss-Prot Acc.P24752.Plays a major role in ketone body metabolism.

Note: This product is for in vitro research use only