

ACSS2 Rabbit mAb

Catalog No: #52120

Package Size: #52120-1 50ul #52120-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	ACSS2 Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal antibody
Clone No.	S01-716
Isotype	Rabbit IgG
Purification	Affinity Purified
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	A synthetic peptide of human ACSS2
Conjugates	Unconjugated
Modification	Unmodification
Other Names	ACSS2; ACSA; ACAS2; AceCS; MYH7B; Acetyl CoA synthetase;
Accession No.	Swiss-Prot:Q9NR19GeneID:55902
Uniprot	Q9NR19
GeneID	55902
Calculated MW	Calculated MW: 79 kDa; Observed MW: 79 kDa
Concentration	0.3 mg/ml
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

WB: 1/1000

Images



Western blot detection of ACSS2 in HeLa cell lysates using ACSS2 Rabbit mAb(1:1000 diluted).Predicted band size:79kDa.Observed band size:79kDa.

Background

Swiss-Prot Acc.Q9NR19.Activates acetate so that it can be used for lipid synthesis or for energy generation.

Note: This product is for in vitro research use only