

## AER61 Rabbit mAb

Catalog No: #52126

Package Size: #52126-1 50ul #52126-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	AER61 Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal antibody
Clone No.	S02-0E5
Isotype	Rabbit IgG
Purification	Affinity Purified
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant protein of human AER61
Conjugates	Unconjugated
Modification	Unmodification
Other Names	Extracellular O-linked N-acetylglucosamine transferase
Accession No.	Swiss-Prot:Q5NDL2GeneID:285203
Uniprot	Q5NDL2
GeneID	285203
Calculated MW	Calculated MW: 62 kDa; Observed MW: 62 kDa
Concentration	0.3 mg/ml
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

## Application Details

WB: 1/1000

## Images



Western blot detection of AER61 in C6,HeLa cell lysates using AER61 Rabbit mAb(1:1000 diluted).Predicted band size:62kDa.Observed band size:62kDa.

## Background

---

Swiss-Prot Acc.Q5NDL2.EGF domain-specific O-linked N-acetylglucosamine transferase

---

Note: This product is for in vitro research use only