

AGXT Rabbit mAb

Catalog No: #52128

Package Size: #52128-1 50ul #52128-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	AGXT Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal antibody
Clone No.	S02-7A4
Isotype	Rabbit IgG
Purification	Affinity Purified
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant protein of human AGXT
Conjugates	Unconjugated
Modification	Unmodification
Other Names	AGT; PH1; SPT; AGT1; SPAT; TLH6; AGXT1
Accession No.	Swiss-Prot:P21549GenelD:189
Uniprot	P21549
GeneID	189
Calculated MW	Calculated MW: 43 kDa; Observed MW: 43 kDa
Concentration	0.3 mg/ml
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

WB: 1/1000-1/5000

Images



Western blot detection of AGXT in mouse liver, rat liver lysates using AGXT antibody. Predicted band size: 43kDa. Observed band size: 43kDa.

Background

Swiss-Prot Acc.P21549.Serine--pyruvate aminotransferase

Note: This product is for in vitro research use only