

ALS2CR8 Rabbit mAb

Catalog No: #52138

Package Size: #52138-1 50ul #52138-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	ALS2CR8 Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal antibody
Clone No.	S04-5G7
Isotype	Rabbit IgG
Purification	Affinity Purified
Applications	WB
Species Reactivity	Human
Immunogen Description	Recombinant protein of human ALS2CR8
Conjugates	Unconjugated
Modification	Unmodification
Other Names	ALS2CR8; NYD-SP24
Accession No.	Swiss-Prot:Q8N187GeneID:79800
Uniprot	Q8N187
GeneID	79800
Calculated MW	Calculated MW: 81 kDa; Observed MW: 81 kDa
Concentration	0.3 mg/ml
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

WB: 1/1000-1/5000;

Images



Western blot detection of ALS2CR8 in 293T lysates using ALS2CR8 antibody. Predicted band size: 81 kDa. Observed band size: 81 kDa.

Background

Swiss-Prot Acc.Q8N187. Acts as a transcriptional activator that mediates the calcium- and neuron-selective induction of BDNF exon III transcription. Binds to the consensus calcium-response element CaRE1 5'-CTATTCGAG-3' sequence.

Note: This product is for in vitro research use only