

## Annexin VII Rabbit mAb

Catalog No: #52146

Package Size: #52146-1 50ul #52146-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	Annexin VII Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal antibody
Clone No.	S07-1D2
Isotype	Rabbit IgG
Purification	Affinity Purified
Applications	WB IHC
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant protein of human Annexin VII
Conjugates	Unconjugated
Modification	Unmodification
Other Names	SNX; ANX7; SYNEXIN
Accession No.	Swiss-Prot:P20073GenelD:310
Uniprot	P20073
GenelD	310
Calculated MW	Calculated MW: 53 kDa; Observed MW: 53 kDa
Concentration	0.3 mg/ml
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


## Application Details

WB: 1/2000; IHC: 1/50;

## Images



Western blot detection of Annexin VII in K562,Rat Brain cell lysates using Annexin VII Rabbit mAb(1:1000 diluted).Predicted band size:53kDa.Observed band size:53kDa.



Immunohistochemistry of Annexin VII in paraffin-embedded Human tonsil using Annexin VII Rabbit mAb at dilution 1/50

## Background

Swiss-Prot Acc.P20073.Calcium/phospholipid-binding protein which promotes membrane fusion and is involved in exocytosis.

Note: This product is for in vitro research use only