

## ARHGAP25 Rabbit mAb

Catalog No: #52157

Package Size: #52157-1 50ul #52157-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	ARHGAP25 Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal antibody
Clone No.	S07-6G6
Isotype	Rabbit IgG
Purification	Affinity Purified
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant protein of human ARHGAP25
Conjugates	Unconjugated
Modification	Unmodification
Other Names	KAIA0053; HEL-S-308
Accession No.	Swiss-Prot:P42331GenelD:9938
Uniprot	P42331
GenelD	9938
Calculated MW	Calculated MW: 73 kDa; Observed MW: 73 kDa
Concentration	0.3 mg/ml
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

## Application Details

WB: 1/1000;

## Images



Western blot detection of ARHGAP25 in Rat brain, Raoms, Jurkat cell lysates using ARHGAP25 Rabbit mAb (1:1000 diluted). Predicted band size: 73kDa. Observed band size: 73kDa.

## Background

---

Swiss-Prot Acc.P42331.GTPase activator for the Rho-type GTPases by converting them to an inactive GDP-bound state.

---

Note: This product is for in vitro research use only