

## Axin 2 Rabbit mAb

Catalog No: #52167

Package Size: #52167-1 50ul #52167-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	Axin 2 Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal antibody
Clone No.	S06-3A8
Isotype	Rabbit IgG
Purification	Affinity Purified
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	A synthetic peptide of human Axin 2
Conjugates	Unconjugated
Modification	Unmodification
Other Names	Axil; Axin like protein; Axin2; Axis inhibition protein 2; Conductin;
Accession No.	Swiss-Prot:Q9Y2T1GeneID:8313
Uniprot	Q9Y2T1
GeneID	8313
Calculated MW	Calculated MW: 94 kDa; Observed MW: 94 kDa
Concentration	0.3 mg/ml
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

## Application Details

WB: 1/1000

## Images



Western blot detection of Axin 2 in K562,C6,3T3,HeLa cell lysates using Axin 2 Rabbit mAb(1:1000 diluted).Predicted band size:94kDa.Observed band size:94kDa.

## Background

---

Swiss-Prot Acc.Q9Y2T1. Inhibitor of the Wnt signaling pathway. Down-regulates beta-catenin. Probably facilitate the phosphorylation of beta-catenin and APC by GSK3B .

---

Note: This product is for in vitro research use only