

BAT5 Rabbit mAb

Catalog No: #52174

Package Size: #52174-1 50ul #52174-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	BAT5 Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal antibody
Clone No.	S01-4D1
Isotype	Rabbit IgG
Purification	Affinity Purified
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	A synthetic peptide of human BAT5
Conjugates	Unconjugated
Modification	Unmodification
Other Names	BAT5; NG26; PP199; D6S82E
Accession No.	Swiss-Prot:O95870GeneID:7920
Uniprot	O95870
GeneID	7920
Calculated MW	Calculated MW: 63 kDa; Observed MW: 55 kDa
Concentration	0.3 mg/ml
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

WB: 1/1000

Images



Western blot detection of BAT5 in HeLa, Jurkat cell lysates using BAT5 Rabbit mAb (1:1000 diluted). Predicted band size: 63 kDa. Observed band size: 55 kDa.

Background

Swiss-Prot Acc.O95870.SitesFeature keyPosition(s)Description

Note: This product is for in vitro research use only