

## BMP2 Rabbit mAb

Catalog No: #52178

Package Size: #52178-1 50ul #52178-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	BMP2 Rabbit mAb
Clone No.	S03-3H9
Isotype	Rabbit IgG
Purification	Affinity Purified
Applications	WB
Species Reactivity	Human,Mouse
Immunogen Description	Recombinant protein of human BMP2
Conjugates	Unconjugated
Modification	Unmodification
Other Names	BDA2; BMP-2A; Bmp2; BMP2A; Bone morphogenetic protein 2A;
Accession No.	Swiss-Prot:P12643GeneID:650
Uniprot	P12643
GeneID	650
Calculated MW	38kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

## Application Details

WB 1:500-1:2000

## Images



Western blot analysis of BMP2 expression in HeLa cell lysate

## Background

Bone morphogenetic proteins (BMPs) are members of the TGF $\beta$  superfamily. BMPs are involved in the induction of cartilage and bone formation. Induces cartilage and bone formation.

---

Note: This product is for in vitro research use only