

CBR3 Rabbit mAb

Catalog No: #52198

Package Size: #52198-1 50ul #52198-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	CBR3 Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal antibody
Clone No.	S06-7B3
Isotype	Rabbit IgG
Purification	Affinity Purified
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	A synthetic peptide of human CBR3
Conjugates	Unconjugated
Modification	Unmodification
Other Names	hCBR3; SDR21C2; HEL-S-25
Accession No.	Swiss-Prot:O75828GeneID:874
Uniprot	O75828
GeneID	874
Calculated MW	Calculated MW: 31 kDa; Observed MW: 31 kDa
Concentration	0.3 mg/ml
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

WB: 1/2000-1/10000

Images



Western blot detection of CBR3 in HeLa lysates using CBR3 antibody. Predicted band size: 31kDa. Observed band size: 31kDa.

Background

Swiss-Prot Acc.O75828.Has low NADPH-dependent oxidoreductase activity towards 4-benzoylpyridine and menadione (in vitro).

Note: This product is for in vitro research use only