

## CD127 Rabbit mAb

Catalog No: #52199

Package Size: #52199-1 50ul #52199-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	CD127 Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal antibody
Clone No.	S08-512
Isotype	Rabbit IgG
Purification	Affinity Purified
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	A synthetic peptide of human CD127
Conjugates	Unconjugated
Modification	Unmodification
Other Names	ILRA; CD127; IL7RA; CDW127; IL-7R-alpha
Accession No.	Swiss-Prot:P16871GenelD:3575
Uniprot	P16871
GenelD	3575
Calculated MW	Calculated MW: 52 kDa; Observed MW: 70 kDa
Concentration	0.3 mg/ml
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

## Application Details

WB: 1/1000-1/5000;

## Images



Western blot detection of CD127 in K562 cell lysates using CD127 Rabbit mAb(1:1000 diluted).Predicted band size:52kDa.Observed band size:70kDa.

## Background

---

Swiss-Prot Acc.P16871.Receptor for interleukin-7. Also acts as a receptor for thymic stromal lymphopoietin (TSLP).

---

Note: This product is for in vitro research use only