

Coilin Rabbit mAb

Catalog No: #52224

Package Size: #52224-1 50ul #52224-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	Coilin Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal antibody
Clone No.	S02-4B9
Isotype	Rabbit IgG
Purification	Affinity Purified
Applications	WB
Species Reactivity	Human
Immunogen Description	A synthetic peptide of human Coilin
Conjugates	Unconjugated
Modification	Unmodification
Other Names	CLN80; p80-coilin
Accession No.	Swiss-Prot:P38432GenelD:8161
Uniprot	P38432
GenelD	8161
Calculated MW	Calculated MW: 63 kDa; Observed MW: 80 kDa
Concentration	0.3 mg/ml
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

WB: 1/1000;

Images



Western blot detection of Coilin in K562 cell lysates using Coilin Rabbit mAb(1:1000 diluted). Predicted band size:63kDa.Observed band size:80kDa.

Background

Swiss-Prot Acc.P38432. Component of nuclear coiled bodies, also known as Cajal bodies or CBs, which are involved in the modification and assembly of nucleoplasmic snRNPs.

Note: This product is for in vitro research use only