

## CYB5R3 Rabbit mAb

Catalog No: #52235

Package Size: #52235-1 50ul #52235-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	CYB5R3 Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal antibody
Clone No.	S01-9F3
Isotype	Rabbit IgG
Purification	Affinity Purified
Applications	WB IHC
Species Reactivity	Human,Rat
Immunogen Description	A synthetic peptide of human CYB5R3
Conjugates	Unconjugated
Modification	Unmodification
Other Names	B5R; DIA1
Accession No.	Swiss-Prot:P00387GenelD:1727
Uniprot	P00387
GenelD	1727
Calculated MW	Calculated MW: 34 kDa; Observed MW: 34 kDa
Concentration	0.3 mg/ml
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


## Application Details

WB: 1/1000-1/5000; IHC: 1/50-1/200

## Images



Western blot detection of CYB5R3 in C6,HeLa cell lysates using CYB5R3 Rabbit mAb(1:1000 diluted).Predicted band size:34kDa.Observed band size:34kDa.



Immunohistochemistry of CYB5R3 in paraffin-embedded Human breast cancer tissue using CYB5R3 Rabbit mAb at dilution 1/50

## Background

Swiss-Prot Acc.P00387.Desaturation and elongation of fatty acids, cholesterol biosynthesis, drug metabolism, and, in erythrocyte, methemoglobin reduction.

Note: This product is for in vitro research use only