

FANCB Rabbit mAb

Catalog No: #52278

Package Size: #52278-1 50ul #52278-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	FANCB Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal antibody
Clone No.	S05-8B7
Isotype	Rabbit IgG
Purification	Affinity Purified
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	A synthetic peptide of human FANCB
Conjugates	Unconjugated
Modification	Unmodification
Other Names	FA2; FAB; FACB; FAAP90; FAAP95
Accession No.	Swiss-Prot:Q8NB91GeneID:2187
Uniprot	Q8NB91
GeneID	2187
Calculated MW	Calculated MW: 98 kDa; Observed MW: 98 kDa
Concentration	0.3 mg/ml
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

WB: 1/1000

Images



Western blot detection of FANCB in 3T3,HeLa cell lysates using FANCB Rabbit mAb(1:1000 diluted).Predicted band size:98kDa.Observed band size:98kDa.

Background

Swiss-Prot Acc.Q8NB91.DNA repair protein required for FANCD2 ubiquitination.

Note: This product is for in vitro research use only