

FLI1 Rabbit mAb

Catalog No: #52284

Package Size: #52284-1 50ul #52284-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

| | |
|-----------------------|--|
| Product Name | FLI1 Rabbit mAb |
| Host Species | Recombinant Rabbit |
| Clonality | Monoclonal antibody |
| Clone No. | S01-0C2 |
| Isotype | Rabbit IgG |
| Purification | Affinity Purified |
| Applications | WB IHC |
| Species Reactivity | Human,Mouse |
| Immunogen Description | A synthetic peptide of human FLI1 |
| Conjugates | Unconjugated |
| Modification | Unmodification |
| Other Names | EWSR2; SIC-1; BDPLT21 |
| Accession No. | Swiss-Prot:Q01543GeneID:2313 |
| Uniprot | Q01543 |
| GeneID | 2313 |
| Calculated MW | Calculated MW: 51 kDa; Observed MW: 51 kDa |
| Concentration | 0.3 mg/ml |
| Formulation | 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA |
| Storage | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles. |


Application Details

WB: 1/1000; IHC: 1/50

Images



Western blot detection of FLI1 in Jurkat,C6 cell lysates using FLI1 Rabbit mAb(1:500 diluted).Predicted band size:51kDa.Observed band size:51kDa.



Immunohistochemistry of FLI1 in paraffin-embedded Human lung cancer tissue using FLI1 Rabbit mAb at dilution 1/50

Background

Swiss-Prot Acc.Q01543.Sequence-specific transcriptional activator (PubMed:24100448, PubMed:26316623, PubMed:28255014). Recognizes the DNA sequence 5'-C[CA]GGAAGT-3'.

Note: This product is for in vitro research use only