

FNDC5 Rabbit mAb

Catalog No: #52285

Package Size: #52285-1 50ul #52285-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	FNDC5 Rabbit mAb
Clone No.	S09-4J9
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	A synthesized peptide derived from human FNDC5
Conjugates	Unconjugated
Modification	Unmodification
Other Names	FNDC 5; FRCP2; irisin;
Accession No.	Swiss-Prot:Q8NAU1GeneID:252995
Uniprot	Q8NAU1
GeneID	252995
Calculated MW	Calculated Mol Weight 24 kDa Observed Mol Weight 24 kDa
Concentration	0.3 mg/ml
Formulation	Rabbit IgG in 10mM phosphate buffered saline , pH 7.4, 150mM sodium chloride, 0.05% BSA, 0.02% sodium azide and 50% glycerol.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

WB 1:1000-1:2000

Images



Western blot analysis of FNDC5 expression in Rat muscle cell lysate.

Background

Contrary to mouse, may not be involved in the beneficial effects of muscular exercise, nor in the induction of browning of human white adipose tissue.

Note: This product is for in vitro research use only