

FUBP1 Rabbit mAb

Catalog No: #52293

Package Size: #52293-1 50ul #52293-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	FUBP1 Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal antibody
Clone No.	S06-4F3
Isotype	Rabbit IgG
Purification	Affinity Purified
Applications	WB IHC
Species Reactivity	Human,Mouse,Rat
Immunogen Description	A synthetic peptide of human FUBP1
Conjugates	Unconjugated
Modification	Unmodification
Other Names	FBP; FUBP; hDH V
Accession No.	Swiss-Prot:Q96AE4GeneID:8880
Uniprot	Q96AE4
GeneID	8880
Calculated MW	Calculated MW: 68 kDa; Observed MW: 74 kDa
Concentration	0.3 mg/ml
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB: 1/1000-1/5000; IHC: 1/50;

Images



Western blot detection of FUBP1/FPB in Jurkat,C6 cell lysates using FUBP1/FPB Rabbit mAb(1:1000 diluted).Predicted band size:68kDa.Observed band size:74kDa.



Immunohistochemistry of FUBP1 in paraffin-embedded Human Cholangiocarcinoma using FUBP1 Rabbit mAb at dilution 1/20

Background

Swiss-Prot Acc.Q96AE4.Regulates MYC expression by binding to a single-stranded far-upstream element (FUSE) upstream of the MYC promoter. May act both as activator and repressor of transcription.

Note: This product is for in vitro research use only