

GCLC Rabbit mAb

Catalog No: #52302

Package Size: #52302-1 50ul #52302-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	GCLC Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal antibody
Clone No.	S03-9D9
Isotype	Rabbit IgG
Purification	Affinity Purified
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant protein of human GCLC
Conjugates	Unconjugated
Modification	Unmodification
Other Names	GCL; GCS; GLCL; GLCLC
Accession No.	Swiss-Prot:P48506GenelD:2729
Uniprot	P48506
GenelD	2729
Calculated MW	Calculated MW: 73 kDa; Observed MW: 73 kDa
Concentration	0.3 mg/ml
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

WB: 1/1000

Images



Western blot detection of GCLC in Rat Brain,C6,3T3,Hela cell lysates using GCLC Rabbit mAb(1:1000 diluted).Predicted band size:73kDa.Observed band size:73kDa.

Background

Swiss-Prot Acc.P48506.Glutamate--cysteine ligase catalytic subunit

Note: This product is for in vitro research use only