

## GLI3 Rabbit mAb

Catalog No: #52307

Package Size: #52307-1 50ul #52307-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	GLI3 Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal antibody
Clone No.	S08-2B9
Isotype	Rabbit IgG
Purification	Affinity Purified
Applications	WB
Species Reactivity	Human
Immunogen Description	Recombinant protein of human Gli3
Conjugates	Unconjugated
Modification	Unmodification
Other Names	GLI3-190;GLI3FL;GLI3-83;GLI3 C-terminally truncated form
Accession No.	Swiss-Prot:Q51S56GenelD:2737
Uniprot	Q51S56
GenelD	2737
Calculated MW	Calculated MW: 170 kDa; Observed MW: 200 kDa
Concentration	0.3 mg/ml
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

## Application Details

WB: 1/1000;

## Images



Western blot detection of Gli3 in 293 cell lysates using Gli3 Rabbit mAb(1:1000 diluted).Predicted band size:170kDa.Observed band size:200kDa.

## Background

---

Swiss-Prot Acc.Q5IS56.Transcriptional activator GLI3

---

Note: This product is for in vitro research use only