

## GOLGA1 Rabbit mAb

Catalog No: #52319

Package Size: #52319-1 50ul #52319-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	GOLGA1 Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal antibody
Clone No.	S07-9C6
Isotype	Rabbit IgG
Purification	Affinity Purified
Applications	WB
Species Reactivity	Human
Immunogen Description	A synthetic peptide of human Golgin 97
Conjugates	Unconjugated
Modification	Unmodification
Other Names	golgin-97
Accession No.	Swiss-Prot:Q92805GeneID:2800
Uniprot	Q92805
GeneID	2800
Calculated MW	Calculated MW: 88 kDa; Observed MW: 88 kDa
Concentration	0.3 mg/ml
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

## Application Details

WB: 1/2000-1/10000

## Images



Western blot detection of Golgin 97 in HeLa cell lysates using Golgin 97 Rabbit mAb(1:1000 diluted). Predicted band size:88kDa.Observed band size:88kDa.

## Background

---

Swiss-Prot Acc.Q92805.Involved in vesicular trafficking at the Golgi apparatus level. Involved in endosome-to-Golgi trafficking.

---

Note: This product is for in vitro research use only