

## GPR37 Rabbit mAb

Catalog No: #52320

Package Size: #52320-1 50ul #52320-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	GPR37 Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal antibody
Clone No.	S02-916
Isotype	Rabbit IgG
Purification	Affinity Purified
Applications	WB
Species Reactivity	Human,Mouse
Immunogen Description	A synthetic peptide of human GPCR GPR37
Conjugates	Unconjugated
Modification	Unmodification
Other Names	PAELR; EDNRBL; hET(B)R-LP
Accession No.	Swiss-Prot:O15354GeneID:2861
Uniprot	O15354
GeneID	2861
Calculated MW	Calculated MW: 67 kDa; Observed MW: 67 kDa
Concentration	0.3 mg/ml
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

## Application Details

WB: 1/1000

## Images



Western blot detection of GPCR GPR37 in Rat Brain lysates using GPCR GPR37 Rabbit mAb(1:1000 diluted). Predicted band size:67kDa.Observed band size:67kDa.

## Background

---

Swiss-Prot Acc.O15354.Receptor for the neuroprotective and glioprotective factor prosaposin. Ligand binding induces endocytosis, followed by an ERK phosphorylation cascade.

---

Note: This product is for in vitro research use only