

GPR182 Rabbit mAb

Catalog No: #52322

Package Size: #52322-1 50ul #52322-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	GPR182 Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal antibody
Clone No.	S02-3E6
Isotype	Rabbit IgG
Purification	Affinity Purified
Applications	WB
Species Reactivity	Human,Mouse
Immunogen Description	Recombinant protein of human GPR182
Conjugates	Unconjugated
Modification	Unmodification
Other Names	AMR; 7TMR; ADMR; AM-R; G10D; L1-R; gamrh; hrhAMR
Accession No.	Swiss-Prot:O15218GeneID:11318
Uniprot	O15218
GeneID	11318
Calculated MW	Calculated MW: 45 kDa; Observed MW: 45 kDa
Concentration	0.3 mg/ml
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

WB: 1/1000

Images



Western blot detection of GPR182 in K562,C6,3T3,HeLa cell lysates using GPR182 Rabbit mAb(1:1000 diluted).Predicted band size:45kDa.Observed band size:45kDa.

Background

Swiss-Prot Acc.O15218.Orphan receptor.

Note: This product is for in vitro research use only