

Growth Hormone Receptor Rabbit mAb

Catalog No: #52323

Package Size: #52323-1 50ul #52323-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	Growth Hormone Receptor Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal antibody
Clone No.	S02-6E6
Isotype	Rabbit IgG
Purification	Affinity Purified
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	A synthetic peptide of human Growth hormone receptor
Conjugates	Unconjugated
Modification	Unmodification
Other Names	GHBP; GHIP
Accession No.	Swiss-Prot:P10912GenelD:2690
Uniprot	P10912
GenelD	2690
Calculated MW	Calculated MW: 72 kDa; Observed MW: 120 kDa
Concentration	0.3 mg/ml
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

WB: 1/1000-1/5000;

Images



Western blot detection of Growth hormone receptor in Rat Brain lysates using Growth hormone receptor Rabbit mAb(1:500 diluted).Predicted band size:72kDa.Observed band size:120kDa.

Background

Swiss-Prot Acc.P10912.Receptor for pituitary gland growth hormone involved in regulating postnatal body growth. On ligand binding, couples to the JAK2/STAT5 pathway .

Note: This product is for in vitro research use only