

HNF 4 alpha Rabbit mAb

Catalog No: #52340

Package Size: #52340-1 50ul #52340-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	HNF 4 alpha Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal antibody
Clone No.	S01-7D7
Isotype	Rabbit IgG
Purification	Affinity Purified
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant protein of human HNF-4-alpha
Conjugates	Unconjugated
Modification	Unmodification
Other Names	Hnf4; HNF-4; MODY1; Nr2a1; Tcf14; TCF-14; Hnf4alpha
Accession No.	Swiss-Prot:P49698GenelD:15378
Uniprot	P49698
GenelD	15378
Calculated MW	Calculated MW: 53 kDa; Observed MW: 53 kDa
Concentration	0.3 mg/ml
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

WB: 1/1000;

Images



Western blot detection of HNF-4-alpha in K562,C6,Hela cell lysates using HNF-4-alpha Rabbit mAb(1:1000 diluted).Predicted band size:53kDa.Observed band size:53kDa.

Background

Swiss-Prot Acc.P49698. Transcriptionally controlled transcription factor. Binds to DNA sites required for the transcription of alpha 1-antitrypsin, apolipoprotein CIII, transthyretin genes and HNF1-alpha. May be essential for development of the liver, kidney and intestine. Miscellaneous Binds fatty acids.

Note: This product is for in vitro research use only