

hnRNP G T Rabbit mAb

Catalog No: #52341

Package Size: #52341-1 50ul #52341-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	hnRNP G T Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal antibody
Clone No.	S05-2H3
Isotype	Rabbit IgG
Purification	Affinity Purified
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant protein of human hnRNP G-T
Conjugates	Unconjugated
Modification	Unmodification
Accession No.	Swiss-Prot:O75526GeneID:27288
Uniprot	O75526
GeneID	27288
Calculated MW	Calculated MW: 43 kDa; Observed MW: 43 kDa
Concentration	0.3 mg/ml
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

WB: 1/1000

Images



Western blot detection of hnRNP G-T in K562,Rat Brain,C6,3T3,Hela cell lysates using hnRNP G-T Rabbit mAb(1:1000 diluted).Predicted band size:43kDa.Observed band size:43kDa.

Background

Swiss-Prot Acc.O75526.RNA-binding motif protein, X-linked-like-2

Note: This product is for in vitro research use only