

## Integrin beta 5 Rabbit mAb

Catalog No: #52373

Package Size: #52373-1 50ul #52373-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	Integrin beta 5 Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal antibody
Clone No.	S04-0G6
Isotype	Rabbit IgG
Purification	Affinity Purified
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant protein of human Integrin beta 5
Conjugates	Unconjugated
Modification	Unmodification
Other Names	Integrin beta-5
Accession No.	Swiss-Prot:P18084GenelD:3693
Uniprot	P18084
GenelD	3693
Calculated MW	Calculated MW: 88 kDa; Observed MW: 88 kDa
Concentration	0.3 mg/ml
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

## Application Details

WB: 1/1000

## Images



Western blot detection of Integrin beta 5 in C6,CHO-K1,HeLa cell lysates using Integrin beta 5 Rabbit mAb(1:500 diluted).Predicted band size:88kDa.Observed band size:88kDa.

## Background

---

Swiss-Prot Acc.P18084.Integrin alpha-V/beta-5 (ITGAV:ITGB5) is a receptor for fibronectin. It recognizes the sequence R-G-D in its ligand.(Microbial infection) Integrin ITGAV:ITGB5 acts as a receptor for adenovirus type C.

---

Note: This product is for in vitro research use only