

IQGAP2 Rabbit mAb

Catalog No: #52374

Package Size: #52374-1 50ul #52374-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	IQGAP2 Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal antibody
Clone No.	S02-117
Isotype	Rabbit IgG
Purification	Affinity Purified
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant protein of human IQGAP2
Conjugates	Unconjugated
Modification	Unmodification
Other Names	Ras GTPase-activating-like protein IQGAP2
Accession No.	Swiss-Prot:Q13576GeneID:10788
Uniprot	Q13576
GeneID	10788
Calculated MW	Calculated MW: 181 kDa; Observed MW: 181 kDa
Concentration	0.3 mg/ml
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

WB: 1/1000;

Images



Western blot detection of IQGAP2 in Ramos lysates using IQGAP2 antibody. Predicted band size: 181kDa. Observed band size: 181kDa.

Background

Swiss-Prot Acc.Q13576.Binds to activated CDC42 and RAC1 but does not seem to stimulate their GTPase activity. Associates with calmodulin.

Note: This product is for in vitro research use only