

LGALS3BP Rabbit mAb

Catalog No: #52397

Package Size: #52397-1 50ul #52397-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	LGALS3BP Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal antibody
Clone No.	S03-5B8
Isotype	Rabbit IgG
Purification	Affinity Purified
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant protein of human LGALS3BP
Conjugates	Unconjugated
Modification	Unmodification
Other Names	90K; M2BP; gp90; CyCAP; BTBD17B; MAC-2-BP; TANGO10B
Accession No.	Swiss-Prot:Q08380GeneID:3959
Uniprot	Q08380
GeneID	3959
Calculated MW	Calculated MW: 65 kDa; Observed MW: 100 kDa
Concentration	0.3 mg/ml
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

WB: 1/1000;

Images



Western blot detection of LGALS3BP in HT1080, HepG2, A549 cell lysates using LGALS3BP Rabbit mAb (1:1000 diluted). Predicted band size: 65 kDa. Observed band size: 100 kDa.

Background

Swiss-Prot Acc.Q08380.Promotes integrin-mediated cell adhesion. May stimulate host defense against viruses and tumor cells.

Note: This product is for in vitro research use only