

## MS4A14 Rabbit mAb

Catalog No: #52434

Package Size: #52434-1 50ul #52434-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	MS4A14 Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal antibody
Clone No.	S02-0B8
Isotype	Rabbit IgG
Purification	Affinity Purified
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant protein of human MS4A14
Conjugates	Unconjugated
Modification	Unmodification
Other Names	Testis development protein NYD-SP21
Accession No.	Swiss-Prot:Q96JA4GeneID:84689
Uniprot	Q96JA4
GeneID	84689
Calculated MW	Calculated MW: 77 kDa; Observed MW: 77 kDa
Concentration	0.3 mg/ml
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

## Application Details

WB: 1/2000-1/10000

## Images



Western blot detection of MS4A14 in Jurkat,Rat Brain,C6 cell lysates using MS4A14 Rabbit mAb(1:1000 diluted).Predicted band size:77kDa.Observed band size:77kDa.

## Background

---

Swiss-Prot Acc.Q96JA4.Membrane-spanning 4-domains subfamily A member 14

---

Note: This product is for in vitro research use only