

MSRB3 Rabbit mAb

Catalog No: #52435

Package Size: #52435-1 50ul #52435-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	MSRB3 Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal antibody
Clone No.	S03-9C6
Isotype	Rabbit IgG
Purification	Affinity Purified
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant protein of human MSRB3
Conjugates	Unconjugated
Modification	Unmodification
Other Names	DFNB74
Accession No.	Swiss-Prot:Q8IXL7GenelD:253827
Uniprot	Q8IXL7
GeneID	253827
Calculated MW	Calculated MW: 21 kDa; Observed MW: 21 kDa
Concentration	0.3 mg/ml
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

WB: 1/1000

Images



Western blot detection of MSRB3 in C6 cell lysates using MSRB3 Rabbit mAb(1:1000 diluted).Predicted band size:21kDa.Observed band size:21kDa.

Background

Swiss-Prot Acc.Q8IXL7.Catalyzes the reduction of free and protein-bound methionine sulfoxide to methionine. Isoform 2 is essential for hearing.

Note: This product is for in vitro research use only