

Myelin Protein Zero Rabbit mAb

Catalog No: #52443

Package Size: #52443-1 50ul #52443-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	Myelin Protein Zero Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal antibody
Clone No.	S02-7A2
Isotype	Rabbit IgG
Purification	Affinity Purified
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant protein of human Myelin Protein Zero
Conjugates	Unconjugated
Modification	Unmodification
Other Names	P0; CHM; DSS; MPP; CMT1; CMT1B; CMT2I; CMT2J; CMT4E; CMTDI3; CMTDID; HMSNIB
Accession No.	Swiss-Prot:P25189GenelD:4359
Uniprot	P25189
GenelD	4359
Calculated MW	Calculated MW: 28 kDa; Observed MW: 28 kDa
Concentration	0.3 mg/ml
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

WB: 1/1000

Images



Western blot detection of Myelin Protein Zero in C6 cell lysates using Myelin Protein Zero Rabbit mAb(1:1000 diluted).Predicted band size:28kDa.Observed band size:28kDa.

Background

Swiss-Prot Acc.P25189.Is an adhesion molecule necessary for normal myelination in the peripheral nervous system. It mediates adhesion between adjacent myelin wraps and ultimately drives myelin compaction.

Note: This product is for in vitro research use only