

## Myosin Ib Rabbit mAb

Catalog No: #52446

Package Size: #52446-1 50ul #52446-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	Myosin Ib Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal antibody
Clone No.	S04-6D4
Isotype	Rabbit IgG
Purification	Affinity Purified
Applications	WB IHC
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant protein of human MYO1B
Conjugates	Unconjugated
Modification	Unmodification
Other Names	MYH-1c;MMI-alpha;MMIa
Accession No.	Swiss-Prot:O43795GeneID:4430
Uniprot	O43795
GeneID	4430
Calculated MW	Calculated MW: 132 kDa; Observed MW: 132 kDa
Concentration	0.3 mg/ml
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


## Application Details

WB: 1/1000-1/5000; IHC: 1/20-1/100

## Images



Western blot detection of MYO1B in C6,3T3,HeLa cell lysates using MYO1B Rabbit mAb(1:1000 diluted).Predicted band size:132kDa.Observed band size:132kDa.



Immunohistochemistry of MYO1B in paraffin-embedded Human lung cancer tissue using MYO1B Rabbit mAb at dilution 1/50

## Background

Swiss-Prot Acc.O43795.Unconventional myosin-Ib

Note: This product is for in vitro research use only